


WANTY™
PTA Balloon Catheter



Application

The PTA Balloon Catheters are indicated for peripheral transluminal angioplasty (PTA) for peripheral vascular system (including the iliac artery, femoral artery, iliac femoral artery, popliteal artery and inferior genicular artery).




Product Specifications (0.014" system)

Technical Specification

- Catheter Design
- Catheter Diameter
- Balloon Coating
- Balloon Marker
- Nominal Pressure
- Compliance
- Recommend Guidewire

Coaxial (OTW)
2.7F - 3.2F (0.87mm-1.08mm)
BSL Hydrophilic
2(zero profile)
6 atm (608.0Kpa)
Semi-compliance
0.014"

- Coaxial design

[illegible]